# STORZ

## KARL STORZ — ENDOSKOPE

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IMAGE1 S™ – TC 200, TC 201, TC 300, TC 301, TC 302, TC 304

### INSTRUCTION MANUAL

IMAGE1 S™ – TC 200, TC 201, TC 300, TC 301, TC 302, TC 304

### MANUAL DE UTILIZARE

IMAGE1 S™ – TC 200, TC 201, TC 300, TC 301, TC 302, TC 304

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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual is intended to aid you in the proper installation, connection, and operation of this device. All required details and all actions on your part are clearly presented and explained. Please read this manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

**i** **NOTE:** This IFU is applicable for IMAGE1 S CONNECT® TC 200 and for IMAGE1 S CONNECT® II TC 201. Only when a special feature of IMAGE1 S CONNECT® TC 200 or IMAGE1 S CONNECT® II TC 201 is described the model number TC 200 or TC 201 is mentioned in the IFU to distinguish between CONNECT and CONNECT II. Otherwise the term 'IMAGE1 S CONNECT®' is used.

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

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



**1.1 Possible combinations for  
IMAGE1 S™**

The IMAGE1 S™ comprises IMAGE1 S CONNECT® (TC 200) or IMAGE1 S CONNECT® II (TC 201) and a minimum of one and maximum of three additional LINK modules. The IMAGE1 S™ H3-LINK (TC 300), IMAGE1 S™ X-LINK (TC 301), IMAGE1 S D3-LINK® module (TC 302) and IMAGE1 S™ 4U-LINK module (TC 304) are available as LINK modules.

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**2 Controls, displays,  
connectors, and their uses**

**2.1 IMAGE1 S CONNECT® TC 200  
in combination with the  
H3-LINK, X-LINK and D3-LINK  
modules (see page IV)**

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output sockets  
(KARL STORZ Communication Bus)
- ⑦ LAN interface
- ⑧ USB interface
- ⑨ HD DVI-D video output sockets
- ⑩ Potential equalization connector
- ⑪ 3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ H3-LINK camera connector
- ⑮ X-LINK camera connector
- ⑯ D3-LINK camera connector

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**2.2 IMAGE1 S CONNECT® II  
TC 201 in combination with the  
4U-LINK, X-LINK and D3-LINK  
modules (see page V)**

- ① Power switch
- ② White balance
- ③ USB port
- ④ Potential equalization connector
- ⑤ LINK inputs
- ⑥ SCB input/output socket  
(KARL STORZ Communication Bus)
- ⑦ LAN interface
- ⑧ USB interface
- ⑨ HD DVI-D video output socket
- ⑩ Potential equalization connector
- ⑪ 12G/3G-SDI video output
- ⑫ Power cord socket
- ⑬ LINK output
- ⑭ 4U-LINK camera connector
- ⑮ X-LINK camera connector
- ⑯ D3-LINK camera connector
- ⑰ DisplayPort video output sockets
- ⑱ EXT SYNC connector (for future use)

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Non-ionizing electromagnetic radiation

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
This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

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<b>1</b>	<b>Images of the equipment</b> .....	<b>IV</b>
1.1	Possible combinations for IMAGE1 S™ .....	IV
<b>2</b>	<b>Controls, displays, connectors, and their uses</b> .....	<b>VI</b>
2.1	IMAGE1 S CONNECT® TC200 in combination with the H3-LINK, X-LINK and D3-LINK modules (see page IV) .....	VI
2.2	IMAGE1 S CONNECT® II TC201 in combination with the 4U-LINK, X-LINK and D3-LINK modules (see page V) .....	VII
<b>3</b>	<b>Symbols employed</b> .....	<b>VIII</b>
3.1	Symbols for operation .....	VIII
3.2	Label and packaging symbols .....	IX
<b>4</b>	<b>General information</b> .....	<b>4</b>
4.1	Description of the devices .....	4
4.2	Property rights .....	4
<b>5</b>	<b>Safety instructions</b> .....	<b>5</b>
5.1	Explanation of warnings and cautions .....	5
5.2	General safety information .....	6
5.3	Indications for use .....	10
5.3.1	Intended use .....	10
5.4	Contraindications .....	10
5.5	Patient profile .....	11
5.6	User profile (specialist physician) .....	11
5.7	User profile (assistant) .....	12
5.8	Safety precautions at the site of installation .....	13
5.9	Safety precautions when operating the device .....	13
<b>6</b>	<b>Initial operation</b> .....	<b>14</b>
6.1	Basic equipment .....	14
6.1.1	TC 200 IMAGE1 S CONNECT® .....	14
6.1.2	TC 201 IMAGE1 S CONNECT® II .....	14
6.1.3	TC 300 IMAGE1 S™ H3-LINK .....	14
6.1.4	TC 301 IMAGE1 S™ X-LINK .....	14
6.1.5	TC 302 IMAGE1 S D3-LINK® .....	14
6.1.6	TC 304 IMAGE1 S™ 4U-LINK .....	14
6.2	Unpacking the equipment .....	14
6.3	Preparing the device for use .....	15
6.3.1	Installation .....	15
6.3.2	Connecting the potential equalization .....	15

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## 4 General information

### 4.1 Description of the devices

The IMAGE1 S™ takes a unique and innovative approach to system conception. The modular design principle allows the user to create a system which is tailored to his needs. One IMAGE1 S CONNECT® or IMAGE1 S CONNECT® II module is at the heart of the IMAGE1 S™ and acts as a communication interface between the image and documentation devices.

For use, one IMAGE1 S CONNECT® or IMAGE1 S CONNECT® II plus at least one LINK module is required which receives and forwards the signal. The Light Source Control and S technologies for image optimization provide users with even more comprehensive support. An intuitive user interface assures simple and quick access as well as rapid understanding of menu navigation.

The system offers numerous state-of-the-art functions. These include:

- IMAGE1 S™ visualization modes to assist the user via contrast enhancement and/or more even brightness perception
- Light Source Control for automatic adjustment of the light source
- Integrated image and video capture for simple documentation

Functions are described in greater detail in § 7.4.

### 4.2 Property rights

Protected by at least one of the following US Patents: 7,821,530; 8,274,559; 7,520,853; 5,716,323; 5,788,688; 5,913,817; 6,397,286; 6,484,221; 7,212,227; US 8,089,509; US 8,199,188. For more information, visit: [www.karlstorz.com](http://www.karlstorz.com).

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**5 Safety instructions****5.1 Explanation of warnings and  
cautions**

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

**WARNING:** A Warning indicates that the personal safety of the patient or physician is at risk.

**CAUTION:** A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

**NOTE:** A Note indicates special information about operating the device or clarifies important information.

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**5.2 General safety information**

**WARNING:** Risk of injury and damage to the products: Failure to observe and follow this instruction manual and all of the instruction manuals of the products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals carefully and always follow the instructions given precisely. Check the functioning of the products used in combination.

**WARNING:** Test this device prior to each surgical procedure to ensure that it functions correctly. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.

**WARNING:** Danger of an electric shock! This device may only be connected to a socket with protective earthing. To reduce the risk of electric shock, keep the device out of reach of patients.

**WARNING:** Grounding reliability can only be achieved when the equipment is connected to a 'Hospital Only' or 'Hospital Grade' outlet (i.e., approved for use in an operating room environment). The plug and cord should be routinely inspected. Do not use if damaged.

**WARNING:** Before each use or after a change of viewing modes/settings, the operator should check to ensure the view observed through the endoscope provides a live image (rather than a stored one) and has the correct image orientation.

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 **WARNING:** Consult and follow the storage, transport and operating conditions prior to storing, transporting or operating, respectively. If the device is stored in a warm, humid environment, ensure that it has suitably equilibrated to the operating environment prior to use.

 **WARNING:** To reduce the risk of electrical shock, do not remove cover of unit. Refer servicing to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

 **WARNING:** Perceived electrical shocks may be experienced when using the camera system in conjunction with electrosurgical units.

 **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

 **WARNING:** Refer to the appropriate section of this manual for validated cleaning and disinfection instructions.

 **CAUTION:** Do not store liquids on or above the unit.

 **CAUTION:** In the USA, federal law restricts this device to sale or use by or on the order of a physician.

 **NOTE:** Disposal of these system products, at the end of their useful life, must be in accordance with local regulations.

 **NOTE:** Do not discard as unsorted municipal waste.

 **NOTE:** Discard as electrical/electronic waste; recycle or reuse accordingly.

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**WARNING:** No modification of this equipment is allowed.

**NOTE:** Consult local authorities for reuse/recycle instructions.

**NOTE:** The presence of maximum mains voltage on SIP/SOPs is not considered a normal condition.

**NOTE:** Additional equipment connected to medical electrical equipment must be able to prove their compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see Clause 16 of the 3rd Ed. of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department. (Standard / directive references: IEC 60601-1+A1+A2: 6.8.2.c, 19.2.b, 19.2.c, IEC 60601-1: 7.9.2.5, 8.1, 16.2.d, MDD 93/42/EEC: Annex I clause 13.6.c).

**NOTE:** For use with the product, KARL STORZ recommends the accessories listed in chapter 11. If third-party accessories are used, the user must check the correct operation himself.

**WARNING:** When performing fluorescence imaging procedures, use suitable endoscopes and recommended light cables.

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-  **WARNING:** To ensure safe operation, do not simultaneously touch the device output connectors and the patient.
-  **WARNING:** When endoscopes are used with energized endoscopically used accessories, the patient leakage currents may be additive. This is particularly relevant when using type CF devices. In this case, a Type CF endoscopically used accessory should be used in order to minimize the total patient leakage current.
-  **WARNING:** During discharge of a defibrillator, the image displayed may take up to 4 seconds to recover and in some cases a power cycle may be required to restore the image.
-  **WARNING:** The device should be placed on a flat surface or in a video cart. Do not block air inlets/outlets as overheating may result in fire hazard or device failure during procedure.
-  **WARNING:** Failure to perform proper white balance of system may result in color distortion.
-  **WARNING:** The operator is responsible to limit physical access to the IMAGE1 S™ from unauthorized persons.
-  **WARNING:** In the case of loss of surgical image for greater than 3 seconds, a power cycle may be required to restore the image.
-  **WARNING:** When using high frequency surgical equipment, keep the working part of the active electrode in the field of view of the operator to avoid accidental high frequency burns. Avoid contact between the active electrode and metal parts of the endoscope.

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**5.3 Indications for use****5.3.1 Intended use**

IMAGE1 S™ is a camera control unit (CCU) for use with camera heads or video endoscopes for the visualization and documentation of endoscopic and microscopic procedures.

**5.4 Contraindications**

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

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**5.5 Patient profile**

Use of the IMAGE1 S™ is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health with regard to the endoscopic intervention and declare it adequate.

**5.6 User profile (specialist physician)**

- Recognized medical skills for endoscopic interventions in the relevant specialization
- Adequate powers of comprehension to assess the situation rationally at all times during the application
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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
**5.7 User profile (assistant)**

- Completed training as medical assistant
- Adequate language skills in at least one of the languages used on the device and in the instruction manual
- Be thoroughly trained in the operation and use of the device
- Knowledge of the contents of the instruction manual
- No physical impairments that could diminish perception of activation and visual message

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### 5.8 Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations. The system is not intended for use in areas where there is a risk of explosion. This means, for example:

Wherever easily combustible and explosive inhalation anesthetics and their mixtures are used, the device may not be operated in the demarcated hazard zone. This also applies to easily combustible and potentially explosive chemicals, e.g., skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for potential equalization. It should be connected in accordance with the applicable national regulations.

### 5.9 Safety precautions when operating the device

It is the user's responsibility to make sure the device is safe and operating properly before using the device.

During treatment with the equipment, the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

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### 6.1 Basic equipment

#### 6.1.1 TC 200 IMAGE1 S CONNECT®

- 1 power cord, length 300 cm
- 1 DVI connecting cable, length 300 cm
- 1 SCB connecting cable, length 100 cm
- 1 3G-SDI Cable  
(only delivered in the USA)
- 1 KARL STORZ USB stick
- 1 Keyboard with touchpad
- 1 instruction manual

#### 6.1.2 TC 201 IMAGE1 S CONNECT® II

- 1 power cord, length 300 cm
- 1 DVI connecting cable, length 300 cm
- 1 SCB connecting cable, length 100 cm
- 1 12G/3G-SDI Cable
- 1 DisplayPort cable, length 300 cm
- 1 KARL STORZ USB stick
- 1 Keyboard with touchpad
- 1 instruction manual

#### 6.1.3 TC 300 IMAGE1 S™ H3-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

#### 6.1.4 TC 301 IMAGE1 S™ X-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

#### 6.1.5 TC 302 IMAGE1 S D3-LINK®

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

#### 6.1.6 TC 304 IMAGE1 S™ 4U-LINK

- 1 power cord, length 300 cm
- 1 LINK cable, 20 cm

### 6.2 Unpacking the equipment

Carefully remove the device and accessories from the packaging. Check for missing items and evidence of damage during shipment. Please file any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these may be useful if the unit has to be transported.

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### 6.3 Preparing the device for use

#### 6.3.1 Installation

- i** **NOTE:** Up to three IMAGE1 S™ LINK modules can be connected to the IMAGE1 S CONNECT®.
- i** **NOTE:** Do not cover the lateral air vents. Maintain a minimum distance of 6 cm to other devices.

Place the IMAGE1 S CONNECT® and the corresponding IMAGE1 S™ LINK modules on a flat surface.

#### 6.3.2 Connecting the potential equalization

- i** **NOTE:** The unit's ground line should be installed by a qualified electrician.

The IMAGE1 S CONNECT® and the IMAGE1 S™ LINK modules are equipped with connectors for potential equalization (see fig. A and B).

1. Connect the potential equalization line to the connectors for potential equalization on all modules.
2. Connect the potential equalization line to the connector in the treatment room.

#### 6.3.3 Connecting the IMAGE1 S™ LINK module with IMAGE1 S CONNECT®

Using the LINK cable supplied, connect the LINK module to the IMAGE1 S CONNECT® (see fig. C).

- i** **NOTE:** The IMAGE1 S CONNECT® inputs are numbered. Please connect the LINK modules in ascending order.

#### 6.3.4 Connecting SCB

Pull back the protection device of the SCB connector and insert the connector into one of the SCB sockets of the TC 200 IMAGE1 S CONNECT® or into the SCB socket of the TC 201 IMAGE1 S CONNECT® II. Connect the other end of the cable to the KARL STORZ-SCB (KARL STORZ Communication Bus) control unit or other SCB units (see KARL STORZ SCB control NEO System Instruction Manual) (see fig. D).

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### 6.3.5 Connecting the monitor

Connect the monitor with the IMAGE1 S CONNECT® or CONNECT II (see fig. E). Depending on the types of connectors provided on the monitor, this connection can be generated using either:

- DVI-D or 3G-SDI output for IMAGE1 S CONNECT® TC 200
- DisplayPort, 12G/3G-SDI or DVI-D output for IMAGE1 S CONNECT® II TC 201

The primary working screen should be connected to DVI1 for TC 200.

The primary working screen should be connected to DP1 for TC 201.

**i NOTE:** The connected monitor must support a resolution of 1920 x 1080 pixels.

**i NOTE:** The connected monitor must support a resolution of 3840 x 2160 pixels to support 4K UHD.

**i NOTE:** A compatible 3D monitor must be connected when using a D3-LINK in conjunction with 3D video endoscopes (see chapter 11.11).

### 6.3.6 Connecting USB devices

The USB ports of the IMAGE1 S CONNECT® are suitable for connecting up the following peripherals (see fig. F):

- Keyboard/Mouse
- Printer
- Footswitch
- USB to ACC adaptor
- USB flash drives, FAT32
- IMAGE1 PILOT TC 014

### 6.3.7 Connecting the power supply

Before connecting the devices to the electrical supply line, check that the values on the power sockets concur with the values on the devices' identification plates.

1. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S CONNECT® (see fig. G).
2. Insert the device plug of the power cord into the power supply cord receptacle of the IMAGE1 S™ LINK modules (see fig. H).

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3. Insert the modules' power cords into the power sockets.

**WARNING:** Always use a hospital grade power cord with this camera.

### 6.3.8 Switching on the IMAGE1 S CONNECT®

Switch on IMAGE1 S CONNECT® (see fig. I).

**NOTE:** The connected LINK modules are turned on automatically.

### 6.3.9 Connecting USB storage media

Connect the USB storage media with the IMAGE1 S CONNECT® to store still images and/or videos (see fig. J).

**NOTE:** The USB storage medium must have a FAT 32 format.

### 6.3.10 Connecting the camera head/ video endoscope

Insert the camera head connector into the receptacle of the relevant LINK module (see fig. K).

**NOTE:** Please verify the values set for the primary and secondary sources as well as the layout beforehand (see p. 23).

**NOTE:** Please note the compatibility list of the respective LINK modules (see chap. 11).

**WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

### 6.3.11 Performing the white balance

**NOTE:** Perform a white balance before every use.

Perform the white balance by pressing the  button on the front of the IMAGE1 S CONNECT® or by pressing the  in the Live Menu.

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**6.3.12 Test for proper functioning**

1. Inspect the IMAGE1 S™ module and the camera head and/or the video endoscope for external damage.
2. Ensure that the camera head cable and/or the video endoscope cable is not damaged or kinked.
3. Point the camera head and/or the video endoscope towards an object and check the quality of the display on the monitor.
4. When using 3D video endoscopes, check to make sure the monitor's display is correct. Please ensure that you use the correct polarized eyewear for the monitor.

**i** **NOTE:** *The device and accessories may only be used once the test for proper functioning has been successfully completed.*

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
## 7 Operating instructions

### 7.1 Basic information on operation

The use of IMAGE1 S™ is divided into three menu stages:

- Dashboard (see 7.2)
- Live menu (see 7.3)
- Setup menu (see 7.4)

#### 7.1.1 Keyboard and mouse

The keyboard and mouse can be used to navigate the menus. A keyboard is required to enter patient data.

#### 7.1.2 Camera head buttons

All IMAGE1 FULL HD camera heads feature three buttons for retrieving programmed functions and for menu control:

⑮ scroll up: Scrolls up the menus, options.

⑰ menu/select: Displays the live menu. Once this menu is displayed, this button activates the highlighted option (identical to ENTER on a computer keyboard).

⑯ scroll down: Scrolls down the menus, options.

### 7.2 Dashboard

The Dashboard displays the devices, data storage devices and LINK modules connected to IMAGE1 S CONNECT® (see fig. L). Unconnected devices are shown in gray.

**NOTE:** The Dashboard is shown when no camera head/video endoscope is connected.

**NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.

**NOTE:** The icons display the current settings.

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**7.2.1 Description of functions of the symbols used**

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
### 7.3 Live Menu

The live menu allows you to choose between different camera functions and change the settings of connected devices (see fig. M). Unconnected devices are shown in gray.

You can access the live menu by pressing the enter button on the keyboard or by pressing the menu button  on the camera head/video endoscope or IMAGE1 PILOT (TC 014) left side button.

-  **NOTE:** The icons in the live menu are displayed as soon as a camera head is connected.
-  **NOTE:** Depending on the configuration of your system, individual icons may be deactivated or unavailable. The configuration can be modified in the Setup Menu.
-  **NOTE:** The icons display the current settings.

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**NOTE:** The S technologies are only available with compatible camera heads and video endoscopes (see chapter 11).

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\* SPECTRA A/B is not available in the USA.

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### 7.4 Setup Menu

In the Setup menu, you can change the general settings of the IMAGE1 S CONNECT® as well as the link-specific settings, and optionally save them as a preset.

You can access the Setup Menu by pressing the setup button on the dashboard as well as by pressing the middle camera head button for 7 seconds in the live menu. Or press Shift-Enter on the keyboard in the setup menu.

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## 7.5 Patient Data Privacy

**i NOTE:** By default the IMAGE1 S™ system does not contain any patient data. The patient data management functions (e.g. enter, access, display, store, or transmit) have to be enabled by the user.

**i NOTE:** By enabling the patient data management functions the user accepts the responsibility for all patient related data managed with the help of IMAGE1 S™ according to applicable national data protection laws.

**i NOTE:** KARL STORZ assumes no legal liability or responsibility for loss, correctness, damage or exposure of patient data managed with the help of IMAGE1 S™.

In the Patient Data Privacy settings, you can configure all settings related to patient information. Patient information of any kind is only available if the corresponding settings have been configured.

The configuration of the Patient Data Privacy Menu is password protected by default. The default password for the settings is: 90290245.

Pressing the Manage button opens the password management submenu (see fig. N).

Here you are offered three options (see fig. O):

- To change the current password
- To reset the password to the default password
- To reset the password protection access

Option 3 is provided to give access to the Patient Data Privacy menu settings if the current password is no longer at hand. **Please note that for security reasons all patient data is removed from the device, and the device is reset to the default PDP settings.**

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
The Increased Patient Data Privacy item provides you with preset cyber security settings. When this setting is enabled (= active), all Patient Data Privacy sub configurations are disabled and set to the default settings. The affected sub configurations and their default values are:

- ▷ Display Duration
  - You can determine how long the the patient data is displayed in the Patient Region of the screen.
  - Default: Never
- ▷ Show Patient Tiles
  - You can determine whether the Patient Tile will be shown on the Dashboard and Live Menu or not.
  - Default: No
- ▷ Enable the printing of Patient Data on Prints
  - You can determine whether prints will be labeled with patient information or not.
  - Default: No
- ▷ Add Patient Names to Files Captured
  - You can determine whether media files stored will be labeled with the patient information.
  - Default: No
- ▷ Patient Data Storage
  - You can specify how long data is stored on the device.
  - Default: 1 day
- ▷ View/Save Logged Events
  - You can view the Audit log by activating the Open Log button. Audit Logs can be exported via the 'Save As...' function.

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The Return for Service item (see fig. P) is not affected by the Increased Patient Data Privacy item – it is always available and functions the same as the reset password protection access. **Upon confirmation, it will remove all Patient Data and reset the device settings to its default patient privacy settings.** This function is recommended for use before returning the device for service.

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**8 Reprocessing**

**8.1 General safety information**

**WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

**CAUTION:** When using cleaning or disinfecting solutions, follow the manufacturer's specifications for exposure time and service life.

**CAUTION:** Do not allow water or any liquid to drip inside the unit(s).

**CAUTION:** Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A list of approved chemicals can be requested from [hygiene@karlstorz.com](mailto:hygiene@karlstorz.com).\*  
\*Not available for U.S. Customers

**CAUTION:** National laws and regulations must be observed.

**NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

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## 8.2 Wipe-down disinfection of devices

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

**i NOTE:** KARL STORZ has qualified the resistance of the surfaces vis-à-vis wipe-down disinfection with the agent Dr. Schumacher – Cleanisept Wipes.

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### 9.1 Maintenance and safety check

#### 9.1.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be obtained from your local representative or from the manufacturer.

#### 9.1.2 Safety check

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, the following checks must be carried out on this device at least once a year by persons who, due to their training, knowledge and practical experience, are able to properly carry out such safety checks and who are not subject to any instructions as far as such checking activities are concerned.

#### Visual inspection:

- Check the device and accessories for any mechanical damage which may impair functionality.
- Check that the labels relevant to safety are legible.

#### Electrical measurements:

- Protective grounding measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Touch Current, Earth Leakage Current, and Patient Leakage currents, measured according to IEC 62353 (current edition), and meets the limits specified in the standard.
- Testing at system level in the cart includes patient leakage testing. The IEC 62353 patient leakage measurement includes testing of special (fault) conditions of mains on Applied Parts, and is conducted with one Type CF applied part measured in turn as per Clause 5.3.4.3.1 & Figure B2 of the IEC 62353, (limit 50 µA).

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When tested according to the IEC 60601-1 3.1 ed. at system level in the cart, it includes total leakage current of testing in normal conditions (limit 50 µA) as per clause # 16.6.3

**Test for proper functioning:**

- A functional check must be carried out in accordance with the instruction manual (see § 6.3.12).
- The safety check must be entered in the unit log-book and the results of the check must be documented.
- If the device does not function reliably and/or safely, it must be repaired.

**9.2 Servicing and repair**

Defective items of equipment must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ genuine KARL STORZ parts only.

**9.3 Disposal**

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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#### 9.4 Repair program

For equipment, individual repair is necessary. To bridge the repair period, you usually receive a unit on loan which you then return as soon as you receive the repaired unit.

In the case of repairs, please contact your local KARL STORZ branch or authorized dealer.

##### **Important information**

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among personnel). We reserve the right to return contaminated instruments or devices to the sender.

Repairs, modifications, or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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### 9.5 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- the electrical installation of the location in which the device is connected and operated conforms to the applicable statutes and standards, and
- the unit has been used in accordance with its operating instructions at all times.

### 9.6 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see 'Subsidiaries' section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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### 10.2 Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

### 10.3 Standard compliance

**According to IEC 60601-1 +A1 +A2, IEC 60601-2-18, UL 60601-1, CSA C22.2 No. 601.1-M90:**

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF defibrillation-proof 
- Type of protection against moisture: IPX 0

**According to IEC 60601-1-2:**

Please observe the information on electromagnetic compatibility in § 12 (p. 59).

### 10.4 Directive compliance for IMAGE1 S™ (TC 200, TC 201, TC 300, TC 301, TC 302 and TC 304)

**According to Medical Device Directive (MDD):** This medical device belongs to Class I

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

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### 10.5 Software ownership and licensing

With respect to KARL STORZ products containing software components, Customer has a non-exclusive, limited, non-transferable license to use the programmed logic, computer programs and/or software supplied by KARL STORZ, in connection with, and incorporated into, KARL STORZ Products System (Software) internally, only in the form in which delivered to Customer and for the sole purpose of operating in accordance with KARL STORZ's written instructions for such products sold to Customer (and for no other product). The Software, and all modifications, enhancements and upgrades thereto, will, at all times, remain the property of KSEA. Customer may not duplicate, copy, reverse-engineer, de-compile, or disassemble the Software or in any way modify the Software. Customer has no right to, and may not, create derivatives of the software, and Customer may not attempt to copy, create or re-create the source code of the Software. Any and all such modifications or enhancements to the Software by Customer, in contravention of this license, will immediately become the sole property of KSEA.

Customer hereby acknowledges and agrees (I) that the purchase of KARL STORZ products does not constitute a sale of the Software, (II) that the Software is the property of KARL STORZ, (III) that Customer neither owns nor acquires any claim or intellectual property right in or to the Software, and (IV) KARL STORZ retains all right, title, and interest in and to the Software, at all times, regardless of the form or media in or on which the original or other copies of the Software may exist. In the event of a failure of Customer, or its agents, employees or representatives, to comply with any terms and conditions of the License herein granted, the License will, without any further action by KARL STORZ or any other party, immediately terminate.

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**Attribution Notice**

This device contains open source software. You can view a list of the open source software included in the product, and the open source software license terms, on the IMAGE1 S CONNECT® device itself. From the Dashboard, select 'Setup/General', then 'Software Licenses'.

To obtain a copy of the source code being made publicly available by KARL STORZ Imaging, Inc. ('KARL STORZ') related to open source software used in this KARL STORZ product, you should send your request, along with product part number and serial number, in writing by email to [opensource@karlstorz.com](mailto:opensource@karlstorz.com) or by postal mail to:

KARL STORZ Imaging, Inc.  
Attention: Open Source Software Inquiries  
1 S. Los Carneros  
Goleta, CA 93117  
USA

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Name	Article number
IMAGE1 S™ H3-Z Three-Chip FULL HD Camera Head	TH 100
IMAGE1 S™ H3-ZI Three-Chip FULL HD Camera Head	TH 101
IMAGE1 S™ H3-Z FI Three-Chip FULL HD Camera Head	TH 102
IMAGE1 S™ H3-P Three-Chip FULL HD Pendelum Camera Head	TH 103
IMAGE1 S™ H3-ZA Three-Chip FULL HD Camera Head	TH 104
IMAGE1 S™ H3-M COVIEW® Three-Chip FULL HD Camera Head	TH 106

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Name	Article number
IMAGE1 H3-Z Three-Chip FULL HD Camera Head	22 2200 55-3
IMAGE1 H3-ZI Three-Chip FULL HD Inline Camera Head	22 2200 53-3
IMAGE1 H3-Z FI Three-Chip ICG FULL HD Camera Head	22 2200 85-3*
IMAGE1 H3-P Three-Chip FULL HD Pendelum Camera Head	22 2200 56-3
IMAGE1 H3-ZA Three-Chip FULL HD Camera Head	22 2200 61-3
IMAGE1 H3-FA Three-Chip FULL HD Camera Head	22 2200 60-3

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**11.3 Compatible camera heads for  
IMAGE1 S™ X-LINK TC 301  
with S technologies**

Name	Article number
IMAGE1 S™ D1, 50 Hz	TH 115
IMAGE1 S™ D1, 60 Hz	TH 116
IMAGE1 S™ HX-FI	TH 112**
IMAGE1 S™ HX-P-FI	TH 113**

**11.4 Compatible camera heads for  
IMAGE1 S™ X-LINK TC 301  
without S technologies**

Name	Article number
IMAGE1 S™ HX	TH 110
IMAGE1 S™ HX-P	TH 111

**11.5 Compatible video endoscopes  
for IMAGE1 S™ X-LINK TC 301  
with S technologies**

Name	Article number (PAL/NTSC)
<b>Specialty BRONCHO</b>	
Video Bronchoscope	11900 BP/BN
Video Bronchoscope HD	11910 T/D
Video Bronchoscope	11910 P/S
<b>Specialty ENT</b>	
Video rhinolaryngoscope	11101 VP/VN
Video rhinolaryngoscope	11101 VPS/VNS
<b>Specialty GASTRO</b>	
Video-Gastroscope 9,3/2,8/1100	13821 PKS/ NKS
Video-Gastroscope Ped./Sten. 5,9/2,0/1100	13820 PKS/ NKS

\* Not available for sale in the USA.

\*\* The camera heads TH 112 and TH 113 do not support the CLARA mode (see page 28).

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**11.6 Compatible video endoscopes  
and camera heads for  
IMAGE1 S D3-LINK® TC 302  
with S technology**

Name	Dia. [mm]	Article number (PAL/NTSC)
TIPCAM®1 S 3D, 0°	10	26605 AA
TIPCAM®1 S 3D, 30°	10	26605 BA
TIPCAM®1 S 3D, 0°	4	7240 AA3D
TIPCAM®1 S 3D, 30°	4	7240 BA3D
TIPCAM®1 S 3D, 0°	4	28164 AA3D
TIPCAM®1 S 3D, 30°	4	28164 BA3D
TIPCAM®1 S 3D, 45°	4	7240 FA3D
TIPCAM®1 S 3D, 45°	4	28164 FA3D
VITOM® 3D		TH 200

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**11.7 Compatible camera heads for  
IMAGE1 S™ 4U-LINK TC 304  
with S technology**

Name	Article number (PAL/NTSC)
IMAGE1 S™ 4U	TH 120

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Bestell-Nr. Order no. Nr. comandă	Artikel	Article	Articol
9826 NB-3D	3D-Monitor, 26"	3D Monitor, 26"	Monitor 3D, 26"
9832 NB-3D	3D-Monitor, 32"	3D Monitor, 32"	Monitor 3D, 32"
TM 323	3D-Monitor, 32"	3D Monitor, 32"	Monitor 3D, 32"
TM 330	3D-Monitor, 32"	3D Monitor, 32"	Monitor 3D, 32"
TM 263	3D-Monitor, 26"	3D Monitor, 26"	Monitor 3D, 26"
TM 350	4K 3D-Monitor, 32"	4K 3D Monitor, 32"	Monitor 4K 3D, 32"

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**12 Information on  
electromagnetic  
compatibility (EMC)**

The IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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System Cables and Maximum Lengths used for EMC Compliance				
Type	Shielded	Max Length [m]	Ferrite	Used for
LINK cable	Yes	2	No	Interconnection of inter-module LINK signals
BNC to BNC	Yes	10	No	Connection of 3G/12G-SDI video signals
DVI	Yes	3	No	Connection of DVI signals
DisplayPort (DP)	Yes	3	No	Connection of DisplayPort signals
4-pin mini-DIN	Yes	1	No	Connection of SCB signals
USB	Yes	2	No	Connection to printer, storage device, keyboard
Power cord	No	3	No	Connection of Device to AC Mains
Ethernet	No	2	No	Connection to LAN (RJ-45) port
IMAGE1 PILOT TC 014	Yes	5	No	Connection to VITOM® 3D
External synchronization cable	No	2.5	No	Connection to compatible light source

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**Information on  
electromagnetic  
compatibility (EMC)**

**WARNING:** The use of an accessory or cable with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the IMAGE1 S™. When using an accessory or cable with IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the IMAGE1 S™, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**NOTE:** The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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**Table 1 – Compliance level for Immunity tests**  
Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The user of IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0 % for 1 cycles @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50 Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V <sub>rms</sub> in ISM bands	3 V <sub>rms</sub> on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 6 V <sub>rms</sub> in ISM bands	

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**Table 3 – Test levels for Radiated and conducted Immunity Tests**  
Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The user of IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup> .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the IMAGE1 S™ is used exceeds the applicable RF compliance level above, the IMAGE1 S™ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the IMAGE1 S™.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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<b>Table 4 – Emission class and group</b> <b>Guidelines and manufacturer's declaration – electromagnetic emissions</b>		
The IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The IMAGE1 S™ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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<b>Table 5</b>			
<b>Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™</b>			
The IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.			
Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.			

**HF surgical instrument compatibility**

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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**13 Appendix**
**13.1 Compatible Veterinary video endoscopes for IMAGE1 S™ X-LINK TC 301 with S technologies**

Name	Article number (PAL/NTSC)
<b>Specialty VET</b>	
Veterinary Video Endoscope	60914 NKS/PKS
Veterinary Video Endoscope	60714 NKS/PKS
Veterinary Video Endoscope	60511 NKS/PKS
Veterinary Video Endoscope	60332 NKS/PKS
Veterinary Video Endoscope	60130 NKS/PKS
Veterinary Video Endoscope	60125 NKS/PKS
Veterinary Video Endoscope	60118 NKS/PKS

**13 Anexă**
**13.1 Videoendoscoape de uz veterinar compatibile cu IMAGE1 S™ X-LINK TC 301 cu tehnologii S**

Nume	Număr articol (PAL/NTSC)
<b>Domeniu de specialitate VET</b>	
Video-endoscop veterinar	60914 NKS/PKS
Video-endoscop veterinar	60714 NKS/PKS
Video-endoscop veterinar	60511 NKS/PKS
Video-endoscop veterinar	60332 NKS/PKS
Video-endoscop veterinar	60130 NKS/PKS
Video-endoscop veterinar	60125 NKS/PKS
Video-endoscop veterinar	60118 NKS/PKS

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postfach 230, 78503 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105  
E-Mail: info@karlstorz.com

KARL STORZ Endoskope Berlin GmbH  
Scharnhorststr. 3  
10115 Berlin/Germany  
Phone: +49 (0)30 30 69090, Fax: +49 (0)30 30 19452

KARL STORZ Endoscopy Canada Ltd.  
7171 Millcreek Drive, Mississauga, Ontario  
L5N 3R3 Canada  
Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.  
2151 East Grand Avenue  
El Segundo, CA 90245-5017, USA  
Phone: +1 424 218-8100, Fax: +1 424 218-8525  
Toll free phone: 800 421-0837 (USA only)  
Toll free fax: 800 321-1304 (USA only)  
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.  
1 South Los Carneros Road  
Goleta, CA 93117, USA  
Phone: +1 805 968-7776, Fax: +1 805 685-2588  
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.  
815 N. W. 57th Avenue, Suite 480  
Miami, FL 33126-2042, USA  
Phone: +1 305 262-8980, Fax: +1 305 262-8986  
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.  
Edificio Atlántic, Oficina 3G  
Calle D e/ 1ra y 3ra  
10400 Vedado, Havana, Cuba  
Phone: +537 836 95 06, Fax: +537 836 97 76  
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.  
Lago Constanza No 326,  
Col. Granada, Del. Miguel Hidalgo  
C.P. 11520 Ciudad de México  
Phone: +52 (55) 1101 1520  
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.  
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,  
CEP-04534-011 São Paulo, Brasil  
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680  
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.  
Zufriategui 627 6° Piso, B1638 CAA – Vicente Lopez  
Provincia de Buenos Aires, Argentina  
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773  
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS  
Stamveien1  
1483 Hagan, Norway  
Phone: +47 6380 5600, Fax: +47 6380 5601  
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB  
Storsåtragränd 14  
127 39 Skårholmen, Sweden  
Phone: +46 8 505 648 00  
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY  
Taivaltie 5  
01610 Vantaa, Finland  
Phone: +358 (0)96824774, Fax: +358 (0)968247755  
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG  
Representation Office  
Kęstućio st. 59 / Lenktoji st. 27  
08124 Vilnius, Lithuania  
Phone: +370 5 272 0448  
Mobile: +370 685 67 000  
E-Mail: info-it-lv@karlstorz.com

KARL STORZ Endoskopi Danmark AS  
Skovlytoften 33  
2840 Holte, Denmark  
Phone: +45 45162600, Fax: +45 45162609  
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.  
415 Perth Avenue, Slough  
Berkshire, SL1 4TQ, United Kingdom  
Phone: +44 1753 503500, Fax: +44 1753 578124  
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B. V.  
Displayweg 2  
3821 BT Amersfoort, Netherlands  
Phone: +31 (0)33 4545890  
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.  
Phone: +31 (0)33 4545890  
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S. A. S.  
12, rue Georges Guynemer, Quartier de l'Europe  
78280 Guyancourt, France  
Phone: +33 1 30484200, Fax: +33 1 30484201  
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH  
Landstraße Hauptstr. 148/1/G1  
1030 Wien, Austria  
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479  
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.  
Parque Empresarial San Fernando  
Edificio Munich – Planta Baja  
28830 Madrid, Spain  
Phone: +34 91 6771051, Fax: +34 91 6772981  
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.  
Via dell'Artigianato, 3  
37135 Verona, Italy  
Phone: +39 045 8222000, Fax: +39 045 8222001  
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.  
Capraška 6  
10000 Zagreb, Croatia  
Phone: +385 1 6406 070, Fax: +385 1 6406 077  
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.  
Cesta v Gorice 34b  
1000 Ljubljana, Slovenia  
Phone: +386 1 620 5880, Fax: +386 1 620 5882  
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.  
ul. Bojkowska 47  
44-100 Gliwice, Poland  
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07  
E-Mail: info-pl@karlstorz.com

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KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & CO. KG  
Interchange 21 Tower, Level 33  
399 Sukhumvit Road,  
North Klongtoey, Wattana,  
10110 Bangkok, Thailand  
Phone: +84 28 3823 8000  
Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPlaza Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd., Hong Kong  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwun Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.  
Beijing Branch,  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jianguomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shanghai Branch  
36F, Tower 1 Grand Gateway  
No. 1 Hongqiao Road, Xuhui District,  
200030, Shanghai, People's Republic of China  
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Chengdu Branch  
Room 805, 8F Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
6100414, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Shenyang Branch,  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.  
Guangzhou Branch  
Room 02B, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China  
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.  
15 Orion Road Lane Cove NSW 2066  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: info@karlstorz.au

www.karlstorz.com

**STORZ**  
KARL STORZ—ENDOSKOPE

**KARL STORZ SE & Co. KG**

Dr.-Karl-Storz-Straße 34  
78532 Tuttlingen

Postfach 230  
78503 Tuttlingen  
Germany

Telefon: +49 (0)7461 708-0  
Telefax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
Web: [www.karlstorz.com](http://www.karlstorz.com)